

simCartRx


SimCart Rx™, built specifically for simulation, provides a realistic experience for medication delivery, along with wireless access to patient records in SimEMR® for bedside charting.

This product is designed to support unit dose medication packaging in a one-drawer-per-patient method and is commonly assigned to a single user for the duration of a training scenario. The system includes a cassette drawer design that allows for quick reloading of supplies such as Demo Dose® from Pocket Nurse®. Our common configuration is six cassette drawers and one full width six inch high supplies drawer. There is an instructor quick load function built into the system that will unlock all of the drawers at once for efficient stocking of supplies.

The electronic locking system provides security for each drawer and is operated by the user via the touch screen or keyboard and mouse. Students will log into the system via a user name and password. Events, such as log in and out, are recorded and scanning of patient medications with the supplied bar code scanner are captured via SimEMR® and the ETC Fusion HD™ recording and debriefing system (optional).

SimCart Rx™ includes the CAPS power™ system that permits full time untethered operation by employing a rechargeable battery pack.

22 annual SimEMR® licenses included (2 instructor, 20 student).



Optional model with barcode printer also available. Please contact your KbPort sales rep for more information.

simCartRx and the power of simplicity Connecting what is important to you...

- Complete drawer system control
- In SimEMR®, users can view medication profiles of simulated patients along with the entire simulated electronic health record
- Instructors can create patient profiles within SimEMR® and subsequently label the SimCart Rx™ drawers with patient names and room numbers
- Students can scan a patient wristband using the barcode scanner available on the SimCart Rx™ which will then automatically open the patient chart within SimEMR®
- Users can navigate to the Medication Administration Record and administer medications such as Demo Dose® via the barcode scanner
- Students can remove the patient's ordered medications from the cassette drawers and scan the medication package barcode. The software will recognize the medication barcode and an order fulfillment screen will automatically populate for the student to then document the administered medication directly within SimEMR®
- Users log in using SimEMR® username and password
- All barcode scanning and drawer access (unlocking/locking) using SimCart Rx™ can be recorded using the ETC Fusion HD™ recording and debriefing system
- Expose students to the use of current technology such as electronic health records and barcode scanning which aim to enhance patient safety
- SimCart Rx™ was designed top to bottom to integrate with SimEMR®, which is specifically designed for educational environments

